

NAMRØL



Vera TP5
Art. n°. TP5

Thank you for choosing NAMROL products.

The equipment you have purchased has been designed and manufactured to meet your present needs and, thanks to its in-built versatility and adaptability, your future needs too.

Before using this unit it is good practice to read this handbook carefully and keep it safely.


This will provide you with a full understanding of the unit, will allow you to make the best use of its potentiality and maintain top-level performance, safety and reliability throughout its working life.

And, whatever circumstance may arise, you can count on the professionalism of your NAMROL dealer, the expertise of his technicians and our reliability as constructors.

NAMROL MEDICAL

TABLE OF CONTENTS

FOOT CONTROL	Pag. 3
MOVEMENTS	Pag. 4
ACCESORIES	Pag. 5
CLEANING AND DISINFECTION	Pag. 6
GENERAL INSTRUCTIONS	Pag. 6
WARNINGS	Pag. 7
WARRANTY	Pag. 7
TECHNICAL SPECIFICATIONS	Pag. 8
SYMBOLE	Pag. 9
IDENTIFICATION	Pag. 9
DIMENSIONS	Pag. 9



Reset



Positioning choice (up)



Positioning choice (down)



For emergency stop press any button.

MOVING HEAD, BODY AND BACK

Actuate side handle and move the bodies to the desired position, release and it is locked.


TABLE MOVEMENTS

To move the table, press up or down key

PO BUTTON


To move the table to initial position, press one time PO.

To stop immediatly, press any button.




BASIC POSITIONS PATIENT


All positions with gas piston movements




POSITION 1: Operate lever 4



POSITION 2: Operate lever 4, 2 and 1



POSITION 3: Operate lever 2 and 1



POSITION 4: Operate lever 2 and knob 3

OPERATE LEVER UP AND DOWN ON BOTH SIDES




ROLL HOLDER SUPPORT

Roll holder for easy placement support with output cap.



TABLE MOVEMENT

The table has a unique fixed swivel wheels displacement to the chassis. To move the table, operate the lever up and the system raises the support legs being wheels on the ground for movement. To lock the movement, operate the lever down.



To clean and disinfect the external surfaces of the table (plastic, painted and upholstered parts) use standard commercially available disinfectants and check that the listed active substances do not exceed maximum concentrations indicated below:

- Glutaraldehyde 2%:10%
- Ethanol 96%:40%
- Formaldehyde:0.01%
- Glyoxal:0.15%
- Propanol:35%

WARNING!

Do not exceed indicated concentrations. Do not use products containing alcohol, ammoniac, abrasive substances or benzene. Wipe disinfected areas with a wet cloth (disinfectants may attack surfaces even where diluted).

Apply products with a soft disinfectant-dampened cloth: do not spray the disinfectant directly onto equipment as it could infiltrate the unit/table and compromise proper operation.

NAMROL. cannot be held liable for any damages deriving from failure to observe the warnings and instructions contained in this handbook.

GENERAL INSTRUCTIONS

- Keep the base constantly clean in order to avoid accumulation of organic deposits within the coverings and consequent stinking odours.
- Absolutely avoid pouring of liquids on the base; in case of accidental pouring, dry as soon as possible.
- Full compliance with
- Medical device Directive 93/42/EEC of 14/06/93 and subsequent amendment 2007/47/EC class 1. ensures maximum safety for the user, the patient and the environment.

Purpose of the table:

The table is intended for professional podiatrists as an aid to in-surgery operations that concern affections of the foot and their treatment. These tables have been designed and constructed for use with patients of all ages.

Neither the table can be used in presence of flammable anaesthetic gas mixtures.

When treating patients with pace-makers or similar devices the operator must take the necessary precautions: refer to the specific literature on the matter.

Do not exceed the loads indicated in the “technical specifications” section.

Do not remove the labels on any unit/table devices. Should these deteriorate replace them.

Power supply:

Make sure that the unit/table electrical power supply is powered via an external differential cut-out switch with a 10A current-carrying capacity and a 10mA trip threshold.

Make sure that table power supply is separate.

Make sure that the table power supply is fitted with an efficient grounding (earth) system and that wiring complies with the standards in force.

Climate:

Under extreme conditions (heat, cold, humidity) it is advisable to let a few hours pass between unpacking the equipment and using it for the first time. This precaution provides the time needed to eliminate any condensation that may have formed inside the packaging.

WARRANTY

The warranty will cover all parts NAMRØL with a manufacturing defect (verified).

- 1.- Damage caused by transport in all its forms.
- 2.- Not cover defects caused by fall, blow or misuse of equipment.
- 3.- Deficiencies in the electrical installation of the clinic or home, such as short circuits or excessive voltage drops, etc.
- 4.- Use outside the conditions specified proper installation and operation in the accompanying user manual.
- 5.- Namrol or its authorized representative will repair or replace at its option, without charge, the components of this unit at its discretion appear to be defective.
- 6.- If you have damaged, altered or removed from the table identifications that this has.
- 7.- This warranty supersedes all other warranties expressed or implied, and expressly so not authorize any other person, partnership or association to take on our own any responsibility for our products. Nor will assume responsibility for any direct or indirect that may occur to the buyer, user or third parties.

Power supply	230 V AC monofase/60Hz *115 VAC (USA/CANADA)
Operating voltage	24 voltios
Lifting motor	6000 N
Gas displacement piston head	300 N
Gas displacement piston central body	800 N
Gas displacement piston leg-rest	800 N
Absortion table	4.0 A
Maximum lifting load.	220 Kg
Table weight	80 Kg
Dimensions	670mm x 2.000mm x 570/930mm H.

SYMBOLE



Alternating Current	
Earth	
Voltage (230 V) *115 VAC (USA/CANADA)	
WARNING:see operating instructions before use	
Medical device Directive 93/42/EEC of 14/06/93 and subsequent amendment 2007/47/EC class 1.	
Customer Support Service (task to be effected only by qualified technical staff).	

The table is identified by the CE, whose position is indicated in the photography. In this plate are given the serial number, year of manufacture, upholstery color and model..

When requesting assistance or ordering spare parts from C.S.S. always quote the equipment serial number.



DIMENSIONS AND WIDTH OF MOVEMENT





NAMRØL

DECLARACIÓN DE CONFORMIDAD UE EU Declaration of Conformity

Reglamento de Productos Sanitarios 2017/745 Medical Devices Regulation 2017/745

Fabricante: NAMROL MEDICAL, S.L.
Manufacturer:
Dirección: C/. Progrés , 24-26
Address: 08850, Gavá (Barcelona) Spain

Número de Registro Único
Single Registration Number (SRN): ES-MF-000006718

DECLARA BAJO SU RESPONSABILIDAD, QUE EL PRODUCTO:
DECLARES UNDER ITS SOLE RESPONSIBILITY THAT THE:

Nombre comercial: <i>Trademark:</i>	VERA				
Producto: <i>Product:</i>	CAMILLA TREATMENT BED				
Código UDI-DI <i>UDI-DI code:</i>	08437023142615	08437023142622	08437023142837	08437023142707	08437023142806
REF:	E2	M2	EP2	TP5	M6
Código UDI-DI <i>UDI-DI code:</i>	08437023142820	08437023142639	08437023142424	08437023142844	
REF:	E3	M3	T3	TP3	
Finalidad prevista <i>Intended purpose:</i>	Examen y tratamiento Examination and treatments				

CUMPLE CON LOS REQUISITOS ESENCIALES DEL REGLAMENTO:
CONFORMS WITH THE ESSENTIAL REQUIREMENTS OF THE REGULATION:

REGLAMENTO (UE) 2017/745
Regulation (EU) 2017/745

Reglamento de Productos Sanitarios
Medical Devices Regulation

Clasificación
Classification

Clase I productos sanitarios
Class I medical devices

Regla
Rule

Regla 1, capítulo III, Anexo VIII
Rule 1, chapter III, Annex VIII

Cumple con las siguientes normas internacionales:

Applied Standards in full or in part:

EN 62366:2009	Dispositivos médicos. Aplicaciones de la ingeniería de aptitud de uso de dispositivos médicos. / <i>Medical devices–Application of usability Engineering to medical devices</i>
ISO 14971:2019	Aplicación de la gestión de riesgos a los productos sanitarios. / <i>Application of risk management to medical devices</i>
UNE-EN ISO 13485:2013	Productos sanitarios. Sistemas de gestión de la calidad. Requisitos para fines reglamentarios. / <i>Medical devices. Quality system control</i>
UNE-EN ISO 15223-1:2013	Productos Sanitarios. Símbolos a utilizar en las etiquetas, el etiquetado y la información a suministrar / <i>Graphical Symbols for use in the labelling of medical devices</i>
UNE EN 1041:2014	Información proporcionada por el fabricante de productos sanitarios / <i>Information supplied by the manufacturer with medical devices</i>
UNE-EN 60601-1 2012–	Parte 1: Requisitos generales para la seguridad básica y funcionamiento esencial. / <i>Part 1: General requirements for basic safety and essential performance</i>
UNE-EN 60601-1-2:2008	Equipos electromédicos. Parte 1-2: Requisitos particulares para la seguridad básica y características de funcionamiento esencial. Norma colateral: Compatibilidad SO 009 SOPORTE CON FRENO electromagnética. Requisitos y ensayos. / <i>Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic compatibility - Requirements and tests</i>

Nombre / *Name*: Guillermo Lorman

Firma / *Signature*:



Fecha / *Date*: 1/02/2024

Gavá, Barcelona

Cargo: Director General

Position: General Manager

NAMRØL

C/ Progrés, 24 - 26
08850 Gavá (Barcelona-Spain).
Tels: 0034 93 637 47 47
e-mail: info@namrol.com
www.namrol.com