988 4805

CYBERTECH

ULTRACLEAN PLUS

INSTRUMENT
DISINFECTION
CONCENTRATE

Instrument disinfection concentrate

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CE 0124

12/0021

EN: Type of application: Concentrate for the manual disinfectant cleaning of general and surgical instruments in an immersion bath procedure. Also suitable for use in an ultrasonic bath (without heat) at room temperature. Not suitable for the final disinfection of semi critical and critical medical devices. Instructions for use: Dilute the concentrate (use dispensing bottles or dispensing equipment) with 10°C to 30°C warm water (at least drinking water quality). In the event of very dirty conditions, first clean the instruments mechanically. Then immerse the opened instruments into the solution. All surfaces and cavities to be disinfected have to be fully the solution. Depending on the degree of contamination, the working solution may need to be replaced (RKI recommendation); at least once per working day. After disinfection, the instruments should be rinsed, ideally with demineralised water, and dried thoroughly. Do not mix with cleaning agents or other disinfectants! The preparation recommendations of the instrument manufacturer are to be observed. Semi-critical medical devices must be prepared using a suitable disinfectant and, where applicable, sterilised for terminal disinfection. Critical medical devices are to be treated with a suitable disinfectant and then sterilised. Very good material compatibility with all common instrument materials (e.g. stainless steel, steel, rubber, diamond-coated drills, different plastics such as PMMA, POM, EDPM). If necessary, a full list of material stability data can be requested. Nevertheless, when using new materials, the compatibility is to be checked in advance. Suitable for use in an ultrasonic bath at room temperature. The working solution can be disposed of in the sewage water system. The dosage quantities and application times should be observed as precisely as possible and not be undercut. Overdoses should also be avoided for reasons of environmental protection and material compatibility. For professional users only. Suitable for dental and medical areas. Application times: Bacteria as per EN 13727, EN 14561 and VAH (Association of Applied Hygiene) (dirty conditions): 1.5% 15min., 1.0% 60min., *Candida albicans* as per EN 13624, EN 14562, VAH (dirty conditions): 1.5% 15 min., 1.0% 60 min., Bacteria and *Candida albicans* in an ultrasonic bath as per EN 14561 and EN 14562 (dirty conditions): 1.0% 5min., "Limited virucidal" (enveloped viruses such as HBV, HIV, HCV) as per DW (German Association for the Control of Virus Diseases) / RKI (Robert Koch Institute), EN 14476 and EN 17111 (dirty conditions): 1.5% 15min., TB and mycobacteria (*M. terrae* and *M. avium*) as per EN 14348 and EN 14563 (clean conditions): 4.0% 30min., 3.0% 60min. Active ingredients: 100 g contain: 4.4 g Didecylmethylammonium chloride, 2.75 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine. Safety instructions: H 302 Harmful if swallowed. H 314 Causes severe skin burns and eye damage. H 410 Very toxic to aquatic life with long lasting effects. Avoid release to the environment. Wear protective gloves/eye protection. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor. Dispose of contents/container in accordance with all local / regional / national / international regulations. Contains didecyl dimethyl ammonium chloride, N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine. 648/2004/CE. Contains 5-15% nonionic surfactants, <5 cationic surfactants, <5% EDTA salts, fragrances. Notes on incidents: All serious incidents related to the device shall be reported by the user/patient to the manufacturer and to the competent authority.

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PRISMAN GmbH