



**SE-709 ERGONOMIC STOOL
TABURETE ERGONÓMICO**

MODEL / MODELO	TYPE / TIPO	DESCRIPTION / DESCRIPCIÓN
ARMIS	STOOL / TABURETE	RODABLE BASIS / BASE RODABLE
STRUCTURE / ESTRUCTURA		UPHOLSTERED / TAPIZADO
CHROMED STEEL / ACERO CROMADO		SKAI-ANTIBACTERIAL / SKAI-ANTIBACTERIAS
LIFTING SYSTEM / ELEVACIÓN		COVERS / CARCASAS
GAS PISTON / PISTÓN A GAS		ABS PLASTIC / PLÁSTICO ABS
WEIGHT / PESO	MAX WEIGHT / CARGA MÁXIMA	BACKREST / RESPALDO
7Kg	120Kg.	ADJUSTABLE / REGULABLE
DIMENSIONS / DIMENSIONES		SEAT HEIGHT / ALTURA ASIENTO
450mm x 450mm x 850mm		h min: 460 - h max:590mm

C/ Progrés 24-26
Tels: 0034 93 811 36 88
e-mail: ingenieria@namrol.com

08850 Gavá (Barcelona-Spain)
Fax: 0034 93 894 87 22
www.namrol.com